


## Create thermometer and link to storage unit

1. From the SIMON homepage, search for and select the Provider/Clinic you would like to use.
2. Next, select the 'Clinic Tools' from the main menu.
3. Select 'Manage Assets' underneath 'Clinic Tools'.



WELCOME TO THE  
**STATEWIDE IMMUNIZATION ONLINE NETWORK**

**Default Provider/Clinic** 1

Provider \*  
ALL CARE FAMILY PRACTICE ▼


Clinic \*  
ALL CARE FAMILY PRACTICE ▼

**News**

[03/29/2019] - Welcome to 19.6 Development

*Welcome to 19.6 Development*  
from [Envision Technology Partners](#)

4. From the 'Manage Assets' page, click on the 'Add Asset' green button in the upper right corner.
5. Select 'Add Thermometer'.



**Manage Assets** ⓘ

4 [Add Asset](#) ▼

Showing 0 to 0 of 0 entries

Name	Type	Status	Action
NO DATA AVAILABLE			

Showing 0 to 0 of 0 entries

← Previous    Next →

**Clinic Notes** [Expand](#) [+ Add](#)

There are currently no notes entered for this clinic

6. Type in the information requested
  - a. The Serial Number/ID field is necessary because when the user uploads temperature readings the Serial Number/ID entered on the screen must be contained in the file name

### Add Thermometer i

---

**Name \***

**Date of Purchase**

**Thermometer Type \***

**Manufacturer \***

**Make \***

**Model \***

**Serial Number/ID \***

**Assigned Storage Unit**

**Calibration Due Every (months)**

**Comments**

7. Assign the thermometer to a storage unit by clicking the 'Assigned Storage Unit' dropdown and selecting the desired location.
8. Click 'Create'.
9. A green success message is displayed confirming your changes were saved.

### Add Thermometer i

---

**Name \***

**Date of Purchase**

**Thermometer Type \***

**Manufacturer \***

**Make \***

**Model \***

**Serial Number/ID \***

**Assigned Storage Unit** 7

**ICE BOX**

**Comments**