


Record Temperature Readings (CTM File Upload)

1. Click on 'Storage Units' from the 'Clinic Tools' dropdown in the main menu.
2. Click on the 'Log' dropdown button.




Storage Units ⓘ

ICE BOX	Status/Type
	Pending Approval/Refrigerator
Next Reading:	Last Reading:
Due Before 10/31/19, 1:38 AM	Submitted On 10/30/19, 1:38 PM
10H	Timeframe 12H
	10/30/19, 1:38 AM - 10/30/19, 1:38 PM
	Excursions 0
	Failures 0

Log ▾

- Log Temp. - Upload CTM
- Log Temp. - Manual
- View Readings
- View Storage Unit

3. Select 'Log Temp. – Upload CTM'.



Upload Temperature Reading ⓘ

Storage Units / ICE BOX (Refrigerator) Thermometer (Type) THERMOMETER (CTM) Unit Status Pending Approval

The associated storage unit is not approved. Vaccines cannot be stored in the unit until approved.

Select the CTM file from your computer and upload it below.

Choose File No file chosen

4. Click 'Choose File' and choose the CTM file.
5. The storage unit's thermometer serial number must be in file name.

Pending Temperature Reading - 08/23/19, 3:19 AM (CTM) Cancel Confirm

Success The Temperature Reading is ready for preview

Storage Units / ICE BOX (Refrigerator) THERMOMETER (CTM) (CTM) Pending

Next Reading: 11/23/19, 12:00 AM - OVERDUE

ALERT: This reading is not available until 'Confirm' is selected. Please review the data below and click 'Confirm' to complete the upload process.

Clinic Notes: There are currently no notes entered for this clinic.

Summary

File Name	CTM Interval	Min/Max Temp
Fridge Temp Log 20191113 1234567890.csv	10M	41.0°F / 41.0°F
Total Timeframe	Excursions	Failures
Total Duration: 11H From: 08/23/19, 4:19 PM To: 08/23/19, 3:19 AM	0	0

Excursions

Showing 0 to 0 of 0 entries

Date Temperature

NO DATA AVAILABLE

Showing 0 to 0 of 0 entries

Inspections

Showing 0 to 0 of 0 entries

Date Temperature

NO DATA AVAILABLE

Showing 0 to 0 of 0 entries

Storage Unit

- Edit Storage Unit

Temperature Readings

- Log Temperature

Associated Thermometer

- Edit Thermometer
- Calibrate Thermometer

- After uploading the file, a Success message will appear.
- Review the temperature reading information and click 'Confirm.'
- If a reading triggers a failure, it will return an error message.
- A success message will appear if submitted correctly.

Success The Temperature Reading has successfully been confirmed and saved.

- After a reading has been submitted, the user has 24 hours to delete it before it stays in the system.

Temperature Reading - 08/23/19, 3:19 AM (CTM) +

Storage Units / ICE BOX (Refrigerator) Thermometer (Type) Unit Status

THERMOMETER (CTM) (CTM) Pending

Summary

File Name	CTM Interval	Min/Max Temp
Fridge Temp Log 20191113 1234567890.csv	10M	41.0°F / 41.0°F
Total Timeframe	Excursions	Failures
Total Duration: 11H From: 08/22/19, 4:19 PM To: 08/23/19, 3:19 AM	0	0

Excursions

Showing 0 to 0 of 0 entries

Date Temperature

NO DATA AVAILABLE

Showing 0 to 0 of 0 entries

Action

- Delete Reading
- New Log Temp. - Upload CTM
- New Log Temp. - Manual

Previous Next